DE	Gebrauchsanweisung	2-3
EN	Instructions for use	4-5
F	Notice d'Utilisation	6-7
TR	Kullanma Talimatı	8-9
NL	Instructies voor gebruik	10-11
ES	Instrucciones de uso	12-13
IT	Istruzioni per l'uso	14-15
PT	Instruções de Uso	16-17
RU	Инструкция по использованию	18-19
UK	Інструкція з використання	20-21
HU	Használati utasítás	22-23
FI	Käyttöohjeet	24-25
EL	Οδηγίες χρήσης	26-27
RO	Instrucțiuni de utilizare	28-29
SK	Návod na použitie	30-31
PL	Informacje dla użytkownika	32-33

Reviscon 1.0% / Reviscon Plus 1.6% / Reviscon Mono 2.0%**Produkt:**

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Reviscon 1.0% / Reviscon Plus 1.6% consists of 2.0 ml of viscoelastic solution in a single-use glass syringe for intraarticular use. Reviscon Mono 2.0% consists of 2.4 ml of viscoelastic solution in a single-use glass syringe for intraarticular use.

Description:

Reviscon 1.0% / Reviscon Plus 1.6% / Reviscon Mono 2.0% is sterile, non-pyrogenic, clear, non-inflammatory, highly purified sodium hyaluronate of high molecular weight, dissolved in a buffered physiological saline solution. The highly-purified sodium hyaluronate is obtained from bacteria by fermentation. Reviscon 1.0% / Reviscon Plus 1.6% / Reviscon Mono 2.0% is a product for the relief of the pain and stiffness of the knee joint and other synovial joints in patients with degenerative and traumatic changes to the synovial joint.

Reviscon 1.0% / Reviscon Plus 1.6% / Reviscon Mono 2.0% contains sodium hyaluronate, sodium chloride, disodium hydrogen phosphate, sodium dihydrogen phosphate and water for injection.

	Reviscon 1.0%	Reviscon Plus 1.6%	Reviscon Mono 2.0%
1.0 mL contains	10.0 mg sodium hyaluronate	16.0 mg sodium hyaluronate	20.0 mg sodium hyaluronate
pH	6.8-7.6	6.8-7.6	6.8-7.6
Osmolality	300-350 mOsm/kg	300-350 mOsm/kg	300-350 mOsm/kg
Molecular weight (App.)	3.0 M Da	3.0 M Da	3.0 M Da

Properties and efficacy:

All synovial joints especially the weight-bearing joints, contain viscoelastic sodium hyaluronate. This substance has lubrication and shock absorbing properties, allowing these joints to move normally and painlessly. In patients with degenerative joint disease (osteoarthritis), the viscoelasticity of the synovial fluid is significantly impaired, causing the mechanical stress on the joint and the breakdown of the articular cartilage to greatly increase resulting in limited and painful joint movement. Intraarticular administration of high purity sodium hyaluronate, which has very good viscoelastic properties, can improve the quality of the joint's lubrication. The lubrication and shock absorbing properties of this product reduce pain and improve joint mobility. This effect may last for several months following recommended treatment cycle of intraarticular injections.

Indications:

Reviscon 1.0% / Reviscon Plus 1.6% / Reviscon Mono 2.0% is indicated for the treatment of pain in osteoarthritis (OA) of the knee in patients who have failed to respond adequately to conservative non pharmacologic therapy and simple analgesics (e.g., acetaminophen (paracetamol)).

Contraindications:

Do not administer to patients with known hypersensitivity to hyaluronate preparations.

Intraarticular injections are contraindicated in cases of present infections or skin diseases in the area of the injection site to reduce the potential for developing septic arthritis.

Side Effects:

Following the use Reviscon 1.0% / Reviscon Plus 1.6% / Reviscon Mono 2.0%, patients may experience local symptoms in the joint being treated (pain, sensation of heat, reddening and swelling). The following adverse events have been reported for similar products: mild or moderate arthralgia, in rare cases skin rash, aseptic joint effusions, pruritus and muscular cramps. Further adverse events that have been observed in very rare cases are: allergic reactions, anaphylactic shock, hemarthrosis, phlebitis, pseudosepsis, severe acute inflammatory reaction (SAIR), nasopharyngitis, joint stiffness, tendonitis, bursitis, fever and myalgia.

Warnings and Precautions:

- Reviscon 1.0% / Reviscon Plus 1.6% / Reviscon Mono 2.0% is intended for single use only. The reuse of the product creates a potential infection risks for patients or users.

- Sodium hyaluronate is manufactured by bacterial fermentation and rigorously purified. However the physician should consider the immunological and other potential risks that can be associated with the injection of any biological material.

- Do not reuse syringe. Any repeat usage of the syringe carries a risk of contamination and infection of the patient.

- Do not re-sterilize the pre-filled syringe. Performance will be impaired. Only the solution is sterile but not the outside of the syringe.

- Do not use if package is damaged or opened.

- Do not use after the expiry date printed on the pack.

- In order to avoid overuse of treated joints patients should be advised to relative rest (but no immobilisation) for 24h after each injection.

- Dispose of the syringe and cannula in accordance with accepted medical practice and applicable national, local and institutional requirements.

- There is no evidence concerning the safety of Reviscon 1.0% / Reviscon Plus 1.6% / Reviscon Mono 2.0% in human pregnancy, lactation and the children under 18 years of age. Administration during pregnancy and lactation is at the discretion of the doctor.

Interaction with other agents:

Sodium hyaluronate is incompatible with quaternary ammonium compounds such as benzalkonium chloride solutions. Therefore, Reviscon 1.0% / Reviscon Plus 1.6% / Reviscon Mono 2.0% pre-filled syringes must never come into contact with surgical instruments rinsed with these solutions.

Dosage and Administration:

FOR INTRAARTICULAR INJECTION. FOR SINGLE USE ONLY.

Reviscon 1.0% / Reviscon Plus 1.6% / Reviscon Mono 2.0% should only be used by a physician for intraarticular injection. Use an appropriate size of the needle (19 to 20 gauge is recommended) and length of the needle depending on the joint to be treated. 3 to 5 injections performed one week apart are recommended (Reviscon 1.0%). 1 to 3 injections performed one week apart are recommended (Reviscon Plus 1.6%). 1-2 injection are recommended (Reviscon Mono 2.0%). Several joints may be treated simultaneously and treatment cycles may be repeated. In order to avoid intraarticular infection strict aseptic injection technique has to be applied. It is recommended that an ice -pack be placed on the joint undergoing treatment for 5 -10 minutes in order to prevent pain and swelling. In the case of effusion accompanied by severe pain the fluid must be removed from the affected joint.

1. Twist and pull out

2. Insert and twist the cannula

WRONG!

(Backstop opening at front)

CORRECT!

(Backstop opening at back)

Storage:

Store Reviscon 1.0% / Reviscon Plus 1.6% / Reviscon Mono 2.0% between 2°C and 25°C. Protect from light and shocks. Do not freeze.

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VSY Biotechnology GmbH