Accutrend® Triglycerides

[REF] 11538144	25 test strips
------------------------------	----------------

English

General information: The determination of triglycerides gives your doctor information about your lipid metabolism. This information serves as a basis for further diagnostic measures, therapeutic decisions or monitoring the progress of the patient. However, self-monitoring is not to be seen as a substitute for regular monitoring by your doctor. The results of self-monitoring should be noted down and discussed with the physician. Test strip for the quantitative determination of triglycerides in fresh or heparinised fresh capillary or heparinised/EDTA venous blood. Use only with the following meters: Accutrend GCT or Accutrend Plus. Suitable for self-testing.

How the test works: Each test strip has a test area containing detection reagents. When blood is applied, a chemical reaction takes place and the test area changes colour. The meter records this change in colour and converts the measurement signal to the displayed result using the data previously entered via the code strip.

Preparing and performing the measurement: In order to avoid erroneous measurements, use only the above specified meters for determinations with Accutrend Triglycerides test strips. Please follow the instructions given in this package insert and in the operating instructions of your meter.

Storage and stability: Do not store the pack at temperatures below 2 °C or above 30 °C. Do not perform triglyceride determinations at temperatures below 18 °C or above 30 °C. If you do, the result obtained may only be used for orientation as it may deviate considerably from the true value depending on the temperature. Always close the container with the original stopper immediately after use to protect the test strips from humidity and direct sunlight. If you use an Accutrend Triglycerides test strip at the same time as an Accutrend Glucose and an Accutrend Cholesterol test strip, make sure that you do not mix up the stoppers. The stoppers on the Accutrend Glucose containers are white; those on the other Accutrend strip vials are blue. If the container stoppers are mixed up, the specified shelf-life can no longer be guaranteed. Never allow liquids or disinfectants to get into the container. The stopper contains a non-toxic silicate-based desiccant. If the desiccant is inadvertently ingested, please drink plenty of water! All components of the pack can be discarded along with household waste or, when used in a laboratory or a doctor's office, according to the specific guidelines. Do not use the test strip after the specified expiration date.

Note the following when collecting blood: Prior to collecting the blood sample, wash your hands with warm water and soap. If you use a hand cream, remember that almost all creams contain fat which can affect the triglyceride reading. If you do use a hand cream, wash your hands very thoroughly and rinse carefully with clean water to obtain a true measurement.

Materials provided: Tube of 25 test strips, 1 code strip.

Materials required (but not provided): One of the above listed meters, Accutrend Control TG1 control solution for performance checks, lancing device (e.g. ACCU-CHEK Safe-T-Pro Plus).

Coding: Whenever a new pack of Accutrend Triglycerides test strips is opened, your meter must be recoded. This procedure is described in the operating instructions of your meter. If the code number on the label of the test strip container and the code number displayed by the meter do not match, the meter will not perform a measurement. Always keep the code strip until all test strips of a pack are used up. Store the code strip outside the test strip container. If the code strip is stored in the container, the test strips may be damaged.

Performance check: The performance of the meter and the test strips and their correct handling must be checked regularly using Accutrend Control TG1 control solution. For further details, please refer to the operating instructions for the meter and the package insert of the Accutrend Control TG1 control solution kit. The acceptable control range for a given test strip lot is given on the enclosed value sheet. The target value corresponds to the mean of the lower and upper value of this control range. Values obtained should fall within the defined limits. Each laboratory must establish corrective measures, if values fall outside the limits. Follow the applicable government regulations and local guidelines for quality control.

Limitations and sources of error: Increased intake of vitamin C when eating fruit and vegetable does not affect the triglyceride detection reaction. The detection reaction can be affected by: taking preparation with large doses of ascorbic acid (vitamin C); intravenous infusion of ascorbic acid (vitamin C); during heparin treatment and total fasting, free glycerol can lead to elevated reading; intravenous administration of lipids (fat emulsions in parenteral nutrition); very high triglyceride values (far above the measuring range of 6.9 mmol/L or 600 mg/dL) may trigger a reading of a value in the high range instead of "HI"- when the values are repeatedly above 5 mmol/L or 450 mg/dL, please schedule an examination by your doctor; hematocrit values higher than 55 %, chlorodiazepoxide, gentisic acid, nitrofurantoin, sulfamethoxazol, trimethoprim.

Measuring range: 0.80-6.86 mmol/L (70-600 mg/dL). Linearity of the method is given in range of 0.80-6.86 mmol/L (70-600 mg/dL).

Normal values: Triglyceride values less than 2.3 mmol/L (200 mg/dL) are considered normal. In some patient groups, the reference ranges may be different; you should discuss this individually with your doctor.

Notes for doctors and laboratories:^{1,2,3}

Applications: For the determination of triglycerides in fresh capillary blood; when disorders of lipid metabolism are suspected; for triglyceride self-testing to monitor therapy and dietary measures; screening. If measurements are performed in several persons, as may be the case in a doctor's office or in screening programmes, the fresh capillary blood may also be applied with heparin-coated capillary pipettes (volume 15 µL to 40 µL).

Test principle: Triglycerides are first cleaved by an esterase into glycerol and free fatty acids. Two additional enzymatic steps convert glycerol into hydroacetone phosphate and hydrogen peroxide. In the presence of peroxidase, the hydrogen peroxide oxidises an indicator to a dye and its concentration is then measured by means of reflectance photometry.

Components per test: CHE (Candida cylindracea) 0.15 U; GK (Bacillus stearothermophilus) 0.36 U; GPO (E. coli rec.) 0.08 U; POD (horseradish) 0.63 U; ATP 12 µg; 4-(4-dimethylaminophenyl)-5-methyl-2-(4-hydroxy-3,5-dimethoxyphenyl)-imidazol-dihydrochloride 27 µg; non reactive components 0.43 mg.

Performance data: The data for Accutrend Triglycerides test strips was determined in a series of tests during evaluation. The majority of the data obtained for the test were within the specified range. *Repeatability (within series imprecision):* CV (coefficient of variation) in the normal range 3.4 %, in the pathological range 3.1 %; sample material: EDTA venous blood.

Reproducibility (day-to-day imprecision): CV 2.3 %, sample material: control solution. *Accuracy (method comparisons; mg/dL; regression equations, n samples, coefficient of correlation r):* y = 0.957 x + 10.42 or 1.052 x – 4.33 (n = 43 -144; r = 0.891). Sample material: fresh capillary blood; comparison method x: GPO-PAP method, Roche Diagnostics, sample material: EDTA venous plasma. *Detection limit (lowest value detected):* 0.80 mmol/L (70 mg/dL).

Last update: 2011–10

Significant additions or changes are indicated by a change bar in the margin.

For literature references, and an explanation of the symbols used, please refer to the end of this method sheet. A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used. For further information, please refer to the appropriate meter manual and the method sheets of all necessary components.

[REF] 11538144	25 tiras reactivas
------------------------------	--------------------

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© 2011 Roche Diagnostics.
ACCUTREND, ACCU-CHEK and COBAS are trademarks of Roche.